

Liftera

Treatment Guide
Video Access Manual



Liftera Treatment Guide

Treatment guide file for access:

Liftera-A-v2_Facial Treatment Guide_210521.pdf




Liftera Treatment Guide




After downloading the Liftera A Treatment Guide, you can access a wide range of recommended treatment methods and protocols per applicator profile and treatment area.




Pen Applicator

Cheek, Jaws, Double-chin Lifting Shot Treatment



How to Treat

Description


Lifting shot technique as line-drawing should be applied to the vector direction instead of rubbing mode. In case of overlapping on the same spot, irradiation starting point should be different with 1-2cm interval to control pain-level and make the thermal energy to be properly spread.


Note







- Do not irradiate below the mouth line (marked), because of nerve location.
- When treating of double-chin, avoid to apply around end of ear area (marked) due to nerve location.
- When treating of double-chin, concentrate on close contact of the cartridge to the skin surface while pressing of the applicator little more than general case.
- Energy from cartridge should not be irradiated in one spot while pressing the Shot button.
- When stopping or ending the procedure, always remove the applicator without pressing the Shot button.
- Adjust the moving speed according to Hz.

Liftera-A 03

Pen Treatment Guide



Treatment Area	Skin Condition	Using cartridge type
Mid-face (Cheek, Chin)	Large face size / having subcutaneous fat, thicker skin type	

	Treatment Area	Cartridge	Energy(J)	Shot number (Both sides)	
STEP 01			0.7~0.9	2,000 Dots Shot	<p>Considering the patient's pain level and the treatment area, irradiate the 2,000 Dots around lateral cheek in the ear direction for lifting by using 4.5 mm Pen applicator. Energy set and its irradiation is recommendable starting from 0.7J and increasing each one level by 0.1 J, considering the patient's pain level. Avoid to treat below mouth line (marked) due to nerve location.</p>
STEP 02			0.7~0.9	2,000 Dots Shot	
STEP 03			0.3~0.4	2,000 Dots Shot	


Total PEN 6,000 shot

Liftera-A 16

Liftera Treatment Guide

However, you may also access video files that provide demonstrations per each treatment method and protocol. Follow these instructions:

1. Download the treatment guide (PDF) and video folders, and they should be located into a single zipped folder.





- chapter01_GuidelinesPerTreatmentAreas
- chapter02_Pen_TreatmentGuide_MidFace
- chapter03_Pen_TreatmentGuide_FullFaceLifting&ToningCare
- chapter04_Pen&Line_TreatmentGuide_MidFace
- chapter05_Pen&Line_TreatmentGuide_FullFaceLifting&ToningCare
- Liftera-A-v2_Facial Treatment Guide_210521.pdf


Liftera Treatment Guide

2. After reviewing the treatment guide located in the folder with the videos, each treatment method and protocol page has a “play video” icon. See below example.

Pen Treatment Guide





Treatment Area	Skin Condition	Using cartridge type
Mid-face (Cheek, Chin)	Large face size / having subcutaneous fat, thicker skin type	



Treatment Area	Cartridge	Energy(J)	Shotnumber (Both sides)







Liftera Treatment Guide

3. Click on the "play video" icon and the below notification will appear on screen. Click "Open" and a default video player will play the selected video demonstration.



Liftera Treatment Guide

4. Be sure to keep the treatment guide PDF file and respective video folders in a single folder, otherwise the play icons and function will not work.

-  chapter01_GuidelinesPerTreatmentAreas
-  chapter02_Pen_TreatmentGuide_MidFace
-  chapter03_Pen_TreatmentGuide_FullFaceLifting&ToningCare
-  chapter04_Pen&Line_TreatmentGuide_MidFace
-  chapter05_Pen&Line_TreatmentGuide_FullFaceLifting&ToningCare
-  Liftera-A-v2_Facial Treatment Guide_210521.pdf