

Zionic
rollactive RF

THE NEW FRONTIER
OF **BODY**
REMODELLING



TERMO SALUD
MEDICAL
& AESTHETIC

zionic
rollactive RF


PRODUCT
PRESENTATION

TERMO SALUD
MEDICAL
& AESTHETIC




TERMO SALUD

Skin ageing



TIGHT SKIN



SKIN WITH SIGNS OF FLACCIDITY




ELASTICITY

CELLULITIS

MUSCLE TONE

HYDRATION FIRMNESS

LOCALISED FAT



→ Treatment zone to **slow down ageing and loss of firmness.**

→ Treatment zone to reduce **unwanted fat accumulation.**


What is ZIONIC? ✓

Monopolar resistive radiofrequency diathermy equipment produced by a handpiece that in turn exerts a rotational massage.


The unique combination of **Deep Rotational Active Massage** and 470 kHz monopolar radiofrequency achieves the most demanding and effective clinical results in the shortest possible time.

MONOPOLAR

It has an emitting electrode and a **neutral return electrode** that closes the electrical circuit.




ERGONOMICS
Easy handling




MARP ENGINE
3 rotation modes
10 speeds
Silent and robust

2 APPLICATORS
Maximum adaptation
to all body areas

DEEP ACTION MAX
ø 19 mm



DEEP ACTION
ø 30 mm



2.8" TOUCH SCREEN
Maximum control of
treatment parameters
Intuitive interface

TOTAL CONTROL
Temperature sensor
real-time monitoring of
treatment temperature

Massage pressure sensor

Safety system for over-
temperature control

TERMO SALUD
MEDICAL
& AESTHETIC



How does it work?

The key to the success of ZIONIC is its **dual thermal and mechanical effect** which guarantees unprecedented clinical efficacy






Effects of MARRP



MARRP
+ RF

DEEP
ACTIVE
ROTATIONAL
MASSAGE

- Destroys adipose tissue
- Activates cell metabolism
- Stimulates collagen and elastin production
- Improves circulation



Indications

Cellulite (grade I, II and III)

Contour reduction

Firming

Body contouring

Muscle toning

Benefits

- Accelerates the production of **collagen and elastin**.
- **Increases the metabolism of fat cells and allows adipose** tissue compaction.
- Promotes the elimination of toxins and waste by **activating circulation**.
- Increases the amount of **nutrients and oxygen to the tissue**.
- Decreases **free radical load**.



Targets

MUSCULAR ACTION



DERMIC ACTION
AND HYPODERMIC



VASCULAR ACTION



EFFECT
TONE-UP




ANTI-CELLULITE



**Body
contouring**

REDUCTION
FROM
ADIPOSIITY



ACTIVATE
CIRCULATION



REAFFIRMING



Cosmetics

Boost ZIONIC results with our **Zionic Triplex Contour Action cosmetics**

A new generation of **post-treatment** cosmetics to enhance and prolong the ZIONIC results even at home.

ZIONIC Triplex Contour Action **combats the problems of skin ageing in the body in a synergistic way.**

With its three lipolytic, draining and firming actions, it combats cellulite, flaccidity and localised fat.




Treatments

Body contouring



Reduces fat deposits
Promotes circulation and
reduces fluid retention
Buttock push-up effect

Firming




Increased elasticity and firmness
Activation of fibroblasts

Anti-cellulite



For cellulite grade I-III
Activation of circulation
Stimulation of fibroblasts
Increased elasticity and firmness

Adiposity reduction



Activates fat metabolism

BEFORE

AFTER



Anti-cellulite treatment

66 year old woman | 8 sessions




Body contouring

38 year old woman : 8 sessions

BEFORE

AFTER



Anti-cellulite treatment

66 year old woman | 8 sessions

Clinical results



ZIONIC has won prestigious international
INTERNATIONAL AWARDS



**IDA
2020**
Industrial design
medical equipment
Bronze



**CLAP
2020**
Industrial design
Miscellaneous
Platinum



**Medical Beauty
Awards 2020**
Radiofrequency
treatment
Silver

Science-based technology

> [Ann Dermatol.](#) 2018 Jun;30(3):296-303. doi: 10.5021/ad.2018.30.3.296. Epub 2018 Apr 23.

Usefulness of Monopolar Thermal Radiofrequency Treatment for Periorbital Wrinkles

Song Hee Han ¹, Yeong Min Yoon ², Yang Won Lee ¹, Yong Beom Choe ¹, Kyu Joong Ahn ¹

Affiliations + expand

PMID: 29853743 PMCID: PMC5929946 DOI: 10.5021/ad.2018.30.3.296

Review > [Semin Cutan Med Surg.](#) 2013 Mar;32(1):9-17.

Noninvasive radio frequency for skin tightening and body contouring

Review > [Int J Hyperthermia.](#) 2008 Feb;24(1):3-15. doi: 10.1080/02656730701769841.

Cellular responses to hyperthermia (40–46 degrees C): cell killing and molecular events

Joseph L Roti Roti ¹

Affiliations + expand


PMID: 18214765 DOI: 10.1080/02656730701769841

> [J Cosmet Dermatol.](#) 2020 Sep;19(9):2317-2324. doi: 10.1111/jocd.13449. Epub 2020 May 13.

Monopolar radiofrequency treatment for facial laxity: Histometric analysis

Dong Hye Suh ¹, Hye-Jin Ahn ², Jong-Kil Seo ², Sang Jun Lee ¹, Min Kyung Shin ², Kye Young Song ³

Professional technology that helps you win customers




Comfortable and painless




Combinable with other therapies to optimise results




Fast and lasting clinical results




Innovative patented technology



High ergonomics, adaptability and handling



Versatility of treatments



Software intuitivo



Technical assistance and continuous training

Partners, KOLs and featured customers




Nuria Soteras



Felicidad Carrera



Carmen Navarro



Dra. Iratxe García

Marketing support program



- Regular online product training
- MULTIMEDIA KIT for social media, product videos and material, before and after pictures
- CLINICAL CASES available
- Organisation of MEETINGS with specialists and opinion leaders


We want you to know...

Termosalud has more than 35 years of experience in the aesthetic sector, which guarantees its quality, professionalism and trust.

Founded in Gijón (Asturias), Termosalud is present **in more than 35 countries** around the world.

A 100% Spanish technology product, developed under the highest quality and safety requirements.

Our qualified team of more than 50 professionals will support you with an agile and constant after-sales service in the development of your activity.



TERMO SALUD

MEDICAL
& AESTHETIC

termosalud.com



certificade
0199.2018



zionic
rollactive RF